Glow Discharge Spectroscopy



GDS900 Glow Discharge Spectroscopy

LECO's Glow Discharge Spectrometer (GDS) offers you state-of-the-art technology designed specifically for routine elemental determination in most conductive ferrous and nonferrous materials. The GDS900 features improved performance, stability, accuracy, and precision in steel, iron (including as-cast), aluminum, copper, zinc, nickel, cobalt, tungsten, and titanium. Low melting alloys, resulfurized steel, powder metals,

LECO

205900

d

and other historically difficult-to-analyze materials can be run in a production environment. Sample types include bulk, sheet, wire, pressed, and mounted samples. Sample holders are available to meet non-standard sample forms.

User-Friendly Cornerstone® Brand Software

LECO's exclusive Cornerstone brand software with touch-screen interface enables complete access to analysis control, method settings, diagnostics, reporting, and more in a highly organized, intuitive, and immersive environment. Designed through a collaboration of customer feedback and innovative engineering, Cornerstone features all of the routine day-to-day operations within a single Analysis screen designed for speed and ease-of-use. Our innovative grouping of sample data into sets and replicates simplifies the data output and automatically calculates relevant statistics, alleviating the need for additional data processing.


0	Cornersto	ne® Gds900			A	ministrator	×	ervice Mode 2/2:	5/2016	10:05 AM	1 ar 👔
0	Ô										
SELCT	MULT	SELECT ALL	907 NOT		ANAEYZ.	MERCE SPLT		UPSE ACL.	REEK		NOE
		C) 2	Nome 854108	Analysis Dat	e Description		Method	Analytes			9
	5	0 2		3:59:01 PM 3/31/2015			Cast Iron NIST 176X and CKD				
	5		BS83H	4.03.46 PM			Cast Iron	2	_	_	_
	2		85181A	3/31/2015 4:08:43 PM			NIST 176X and CKD Gest Iron				
۲ Ç	2	2	MBH 13XPH3-G	3/31/2015 4.13.34 PM			NIST 176X and CKD Cast Iron	De De			
• Q	3)	0	Registry	6/18/2015 11:00:32 AM	Data collec Coefficient	tion to set Registry s	Instrument Conditioning				
• Ç)	0	Warmup	6/18/2015 1.53.18 PM			CKD 24x Single Wavelength	29712731 %	3	M.	
2)-	CR 3	CKD/241	6/18/2015 3:43:59 PM			CKD 24x Single Wavelength	0.0164 % %	168 %	3.23 %	0.00439 %
•)-		СК0/242	6/18/2015 3.46.05 PM			CKD 24x Single Wavelength	0.0815 % 0.00887	s 2.22 %	2.87 %	0.0366 %
1)-	<u>(</u> 8 3	CKD/243A	6/18/2015 3:48:09 PM			CKD 24x Single Wavelength	0.0166 % 0.00790	1. K 2.37 K	2.40 %	0.178 %
· C)-	c8 🕤	CKD244	6/18/2015 3.50.13 PM			CKD 24x Single Wavelength	0.0537 % 0.0839 %	E 2.58 %	2.10 %	0.0268 %
C	5-	CR 3	CKD/245A	6/18/2015 3:52:15 PM			CKD 24x Single Wavelength	0.0185 % 0.00529	3	1.57 %	0.407 %
. 6	5	<u>C</u> R 3	CKD/246	6/18/2015 3.54:19 PM			GKD 24x Single	0.00720 % %	16	0.663 %	0.607 %
2	5	c8 3	14440	6/18/2015 3:56:22 PM			Wavelength CKD 24x Single Wavelength	N E	jie –	-	
2	5	<u></u>	CKD/247	3:56:22 PM 6/18/2015 3:58:25 PM			CKD 24x Single	0.0110 % 3000123	<u>а</u>	1.14 %	0.103 %
2	<		CKD/252				Wavelength CKD 24x Single	-0.177 % 0.0320 %	3	32	-0.00645 %
1	2	0 5	Warmup	7/9/2015 2:37:42 PM			Wavelength	0.0146 % %	0.212 %	0.894 %	0.0445 %
LOGI	N	REPLICATES			+ CA +	SAMPLE		METHO	0	Fe_1000_	30
OFF	0	Average	Std. Dev.	RSD	1	2	3	4	5	5	
	AI	0 0293 %	0.0276 %	94.4	0.0164 %	0 0815 %	0.0166 %	0 0537 %	0.018		18 A.
-,-	в	0.0145 %	0.0310 %	213	0.000585 %	0.00887 %	0.00790 %	0.0839 %	0.005	29 %	-
				0.0455			-				
					1	ê .	3	-	é	7	
				-0.0164							
×.*		2.53 %	0.484 %	0.0164	1.68 %	2.22 %	2.37 %	2.58 %	2.93	1%	
میت شید		2.53 % 2.00 %	0.484 %		1.68 % 3.23 %	2.22 % 2.87 %	2.37 % 2.40 %	2.58 % 2.10 %			
									2.93		
		2.00 %	0.927 %	19.1 46.4	3.23 %	2.87 %	2 40 %	2.10 %	2.93 1.57	'% 7%	
		2.00 % 0.195 %	0.927 % 0.228 %	19.1 46.4 117	3.23 % 0.00439 %	2.87 % 0.0366 %	2 40 % 0.178 %	2_10 % 0.0268 %	2.93 1.57 0.40	'% 7% 19%	
1445) 1447 1447 1447	Si P S Mo	2 00 % 0.195 % 0.0481 % 0.322 %	0 927 % 0.228 % 0.0460 % 0.434 %	19.1 46.4 117 95.7 135	3.23 % 0.00439 % 0.136 % 0.577 %	2.87 % 0.0366 % 0.0477 % 1.20 %	2 40 % 0.178 % 0.0793 % 0.277 %	2.10 % 0.0268 % 0.0182 % 0.0710 %	2.93 1.57 0.40 0.034	'% 7% 19%	
	Si P S Mo Ti	2.00 % 0.195 % 0.0481 % 0.322 % 0.0742 %	0.927 % 0.228 % 0.0460 % 0.434 % 0.110 %	19.1 46.4 117 95.7 135 148	3 23 % 0.00439 % 0 136 % 0.577 % 0.00166 %	2.87 % 0.0366 % 0.0477 % 1.20 % 0.314 %	2.40 % 0.178 % 0.0793 % 0.277 % 0.0233 %	2.10 % 0.0268 % 0.0182 % 0.0710 % 0.0324 %	2.93 1.57 0.40 0.034 0.095 0.007	* % 7 % 19 % 10 % 70 %	
ngti Ngti Ngti Ngti	Si P S Mo Ti As	2 00 % 0.195 % 0.0481 % 0.322 % 0.0742 % 0.0205 %	0 927 % 0 228 % 0 0460 % 0 434 % 0 110 % 0 00779 %	19.1 46.4 117 95.7 135 148 38.0	3 23 % 0.00439 % 0.136 % 0.577 % 0.00166 % 0.0206 %	2 87 % 0.0366 % 0.0477 % 1.20 % 0.314 % 0.0299 %	2.40 % 0.178 % 0.0793 % 0.277 % 0.0233 % 0.0308 %	2 10 % 0.0258 % 0.0182 % 0.0710 % 0.0324 % 0.0219 %	2.93 1.57 0.40 0.034 0.095 0.067 0.016	*% 7% 19% 0% 3%	
1445) 1447 1447 1447	Si P S Mo Ti As Sn	2.00 % 0.195 % 0.0481 % 0.322 % 0.0742 % 0.0205 % 0.0601 %	0 927 % 0 228 % 0 0460 % 0.434 % 0.110 % 0 00779 % 0.0606 %	19.1 464 117 957 135 148 380 101	3 23 % 0.00439 % 0.136 % 0.577 % 0.00166 % 0.0206 % 0.0206 %	2 87 % 0.0366 % 0.0477 % 1.20 % 0.314 % 0.0314 % 0.0654 %	2.40 % 0.178 % 0.0733 % 0.277 % 0.0233 % 0.0308 % 0.116 %	2 10 % 0.0268 % 0.0182 % 0.0710 % 0.0324 % 0.0329 % 0.161 %	2.90 1.57 0.402 0.095 0.087 0.016 0.057	* % 7 % 19 % 10 % 3 % 6 %	
ngti Ngti Ngti Ngti	Si P S Mo Ti As	2 00 % 0.195 % 0.0481 % 0.322 % 0.0742 % 0.0205 %	0 927 % 0 228 % 0 0460 % 0 434 % 0 110 % 0 00779 %	19.1 46.4 117 95.7 135 148 38.0	3 23 % 0.00439 % 0.136 % 0.577 % 0.00166 % 0.0206 %	2 87 % 0.0366 % 0.0477 % 1.20 % 0.314 % 0.0299 %	2.40 % 0.178 % 0.0793 % 0.277 % 0.0233 % 0.0308 %	2 10 % 0.0258 % 0.0182 % 0.0710 % 0.0324 % 0.0219 %	2.93 1.57 0.40 0.034 0.095 0.067 0.016	* % 7 % 19 % 10 % 3 % 6 %	
ngti Ng ^a Ng ^a	Si P S Mo Ti As Sn Pb	2.00 % 0.195 % 0.0481 % 0.322 % 0.0742 % 0.0205 % 0.0601 %	0 927 % 0 228 % 0 0460 % 0.434 % 0.110 % 0 00779 % 0.0606 %	19.1 46.4 117 95.7 135 148 38.0 101 32.6	3 23 % 0.00439 % 0.136 % 0.577 % 0.00166 % 0.0206 % 0.0206 %	2.87 % 0.0366 % 0.0477 % 1.20 % 0.314 % 0.0299 % 0.0634 % 0.0136 %	2 40 % 0.178 % 0.0733 % 0.277 % 0.0233 % 0.0308 % 0.116 % 0.00986 %	2 10 % 0.0268 % 0.0162 % 0.0710 % 0.0324 % 0.0219 % 0.161 % 0.0103 %	2.93 1.57 0.40 0.095 0.057 0.016 0.057 0.0069	* % 7 % 19 % 10 % 70 % 70 % 76 % 54 %	
ngti Ng ^a Ng ^a	Si P S Mo Ti As Sn Pb	2 00 % 0.195 % 0.0481 % 0.322 % 0.0742 % 0.0742 % 0.0206 % 0.0601 %	0 927 % 0 228 % 0 0460 % 0.434 % 0.110 % 0 00779 % 0.0606 %	19.1 46.4 117 95.7 135 148 38.0 101 32.6	3 23 % 0.00439 % 0.136 % 0.577 % 0.00166 % 0.0206 % 0.0124 % 0.00124 %	2.87 % 0.0366 % 0.0477 % 1.20 % 0.314 % 0.0299 % 0.0634 % 0.0136 %	2.40 % 0.178 % 0.0793 % 0.0277 % 0.0233 % 0.0308 % 0.0308 % 0.116 % 0.00906 %	2 10 % 0.0268 % 0.0162 % 0.0710 % 0.0324 % 0.0219 % 0.161 % 0.0103 %	2.93 1.57 0.402 0.034 0.095 0.067 0.066 0.057 0.0063	* % 7 % 19 % 10 % 70 % 54 %	
ayt Ayt Ayt Ayt	Si P S Mo Ti As Sn Pb	2 00 % 0.195 % 0.0481 % 0.322 % 0.0742 % 0.0742 % 0.0206 % 0.0601 %	0 927 % 0 228 % 0 0460 % 0.434 % 0.110 % 0 00779 % 0.0606 %	19.1 46.4 117 95.7 135 148 38.0 101 32.6	3 23 % 0.00439 % 0.136 % 0.577 % 0.00166 % 0.0206 % 0.0124 % 0.00124 %	2.87 % 0.0366 % 0.0477 % 1.20 % 0.314 % 0.0299 % 0.0634 % 0.0136 %	2.40 % 0.178 % 0.0793 % 0.0233 % 0.0233 % 0.0306 % 0.116 % 0.00996 %	2 10 % 0.0268 % 0.0162 % 0.0710 % 0.0324 % 0.0219 % 0.161 % 0.0103 %	2.93 1.57 0.402 0.034 0.095 0.067 0.066 0.057 0.0063	* % 7 % 19 % 0 % 3 % 6 % 54 %	
	Si P S Mo Ti As Sn Pb	2.00 % 0.195 % 0.0481 % 0.0322 % 0.0205 % 0.0205 % 0.00082 % VOIDCE	0 927 % 0 228 % 0 0460 % 0.434 % 0.110 % 0 00779 % 0.0606 %	19.1 46.4 117 95.7 135 148 38.0 101 32.6	3 23 % 0.00439 % 0.136 % 0.577 % 0.00166 % 0.0206 % 0.0124 % 0.00124 %	2.87 % 0.0366 % 0.0477 % 1.20 % 0.314 % 0.0299 % 0.0634 % 0.0136 %	2.40 % 0.178 % 0.0793 % 0.0233 % 0.0233 % 0.0308 % 0.116 % 0.00906 %	2 10 % 0.0268 % 0.0162 % 0.0710 % 0.0324 % 0.0219 % 0.161 % 0.0103 %	2.93 1.57 0.402 0.034 0.095 0.067 0.066 0.057 0.0063	* % 7 % 19 % 0 % 3 % 6 % 54 %	
	Si P S Mo Ti As Sn Pb	2 00 % 0.195 % 0.0481 % 0.322 % 0.0742 % 0.0742 % 0.0206 % 0.0601 %	0.927 % 0.228 % 0.0460 % 0.454 % 0.454 % 0.00779 % 0.0606 % 0.00260 %	19.1 46.4 117 95.7 135 148 38.0 101 32.6	3 23 % 0.00439 % 0.136 % 0.577 % 0.00166 % 0.0206 % 0.0124 % 0.00124 %	2.87 % 0.0366 % 0.0477 % 1.20 % 0.314 % 0.0299 % 0.0634 % 0.0136 %	2.40 % 0.178 % 0.0793 % 0.0233 % 0.0233 % 0.0308 % 0.116 % 0.00906 %	2.10 % 0.0258 % 0.0102 % 0.0710 % 0.0324 % 0.0219 % 0.161 % 0.0103 %	2.93 1.57 0.400 0.034 0.057 0.016 0.057 0.016 0.057 0.016 0.057 0.016 0.057	* % 7 % 19 % 0 % 3 % 6 % 54 %	
	Si P S Mo Ti As Sn Pb	2.00 % 0.195 % 0.0481 % 0.0322 % 0.0205 % 0.0205 % 0.00082 % VOIDCE	0.927 % 0.228 % 0.0460 % 0.454 % 0.454 % 0.00779 % 0.0606 % 0.00260 %	19.1 46.4 117 95.7 135 148 38.0 101 32.6 ECC.	3 23 % 0.00439 % 0.136 % 0.577 % 0.00166 % 0.0206 % 0.0124 % 0.00124 %	2.87 % 0.0366 % 0.0477 % 1.20 % 0.314 % 0.0299 % 0.0634 % 0.0136 %	2.40 % 0.178 % 0.0793 % 0.0233 % 0.0306 % 0.116 % 0.00996 %	2.10 % 0.0258 % 0.0102 % 0.0710 % 0.0324 % 0.0219 % 0.161 % 0.0103 %	2.93 1.57 0.400 0.034 0.057 0.016 0.057 0.016 0.057 0.016 0.057 0.016 0.057	* % 7 % 19 % 7	
	Si P S Mo Ti As Sn Pb	2.00 % 0.195 % 0.0481 % 0.0322 % 0.0205 % 0.0205 % 0.00082 % VOIDCE	0.927 % 0.228 % 0.0460 % 0.454 % 0.454 % 0.00779 % 0.0606 % 0.00260 %	19.1 46.4 117 95.7 135 140 38.0 101 32.6 0C.C. 0C.C.	3 23 % 0.00439 % 0.136 % 0.577 % 0.00166 % 0.0206 % 0.0124 % 0.00124 %	2.87 % 0.0366 % 0.0477 % 1.20 % 0.314 % 0.0299 % 0.0634 % 0.0136 %	2.40 % 0.178 % 0.0793 % 0.0233 % 0.0306 % 0.116 % 0.00996 %	2.10 % 0.0258 % 0.0102 % 0.0710 % 0.0324 % 0.0219 % 0.161 % 0.0103 %	2.93 1.57 0.400 0.034 0.057 0.016 0.057 0.016 0.057 0.016 0.057 0.016 0.057	* % 7 % 19 % 7	

Software Features & Benefits


The software is divided into four main sections–Analysis, Diagnostics, Settings, and Instrument–for simplified navigation and organization. Toolbars, sliders, and drop-down menus make it easy to set parameters for calibration and data processing. The software also included real-time monitoring of ambient parameters and spectra, with a fully animated system diagram.

Advanced interactive diagnostic features include a thorough digital on-board manual, maintenance animations, photo illustrations, and screen captures that quickly provide the direction needed without having to refer to multiple manuals.

Cornerstone also supports a multilingual interface, user permissions, extended data archiving and filtering, compatibility with various Laboratory Information Management Systems (LIMS), and flexible reporting capabilities.



Instrument Highlights and Features



Glow Discharge Source Advantages

- Simple, Linear Calibrations: GD-AES provides narrow emission lines, less interference, increased dynamic range, and less complex spectra when compared to other sources.
- Controlled Excitation: GD-AES provides a non-thermal source to uniformly remove sample material from the surface. Excitation of the sample occurs away from the surface, reducing chemical and metallurgical effects. Very little sample-to-sample carryover allows quick matrix changes. Uniform sample excitation offers improved precision.
- Interchangeable Anodes: Choice of 4 mm and 2 mm anodes to optimize spot size to sample size and type.
- Reduced reference material consumption allows more acquisition between required resurfacing and shallower burn spots requiring less material removal during resurfacing.

Detection System Advantages

- The inherent sensitivity, dynamic range, and linearity of CCD detectors coupled with the GD source ensures stability, flexibility, and performance for all bulk elemental applications.
- The all CCD detection system provides full wavelength coverage from 160 nm to 460 nm.
- 50 pm (0.050 nm) resolution to differentiate even the most complex features of bulk elemental spectra.

	1800
	GDS900
U	J



System Level User Advantages

- Very easy to operate, maintain and service. Minimal training required.
- Low consumable consumption for minimized cost per test.
- Automatic cleaning between samples saves time, minimizes matrix effects for increased precision.
- Open lamp design and ambidextrous controls allow both right and left hand operation.

Reliability with Trusted Service and Support

- Knowledgeable sales force with a customer-centered focus dedicated to helping you understand and identify the best instrumentation fit for your application.
- State-of-the-art Technical Services Laboratory with experienced technical application chemists to assist in method development and other application related requests.
- Global and regional LECO service network comprised of regional support centers and over 25 international LECO subsidiaries, dedicated to providing service and support offerings, including field service visits over the lifetime of the instrument.



Model Availability

GDS900C GDS900 Package with PC, touch-screen, software

GDS900RC GDS900 Package with, PC, touch-screen, software, nitrogen purge recirculation

Options



Enhance the ergonomics of the system with the choice of an optional integrated desk or mobile workstation.



Our optional extension spectrometer extends spectral coverage from 460 to 850 nm, allowing access to higher wavelength elements like sodium, potassium, and lithium.



Optional recirculating nitrogen purge system reduces nitrogen consumption by ~1200% annually.



Optional vacuum pump sound abatement enclosure reduces ambient noise from vacuum pump by $\sim 7 \text{ dBa}$.

LECO—Your source for total analytical solutions



CS Series: Carbon/Sulfur by Combustion

- Carbon and sulfur determination by combustion infrared detection
- Quick, accurate, and affordable determination for both production control and research
- Calibration, analysis, evaluation, and diagnostic functions accessible via user-friendly Cornerstone brand software

ONH Series: Oxygen/Nitrogen/Hydrogen by Fusion

- Oxygen, nitrogen, and hydrogen by inert gas fusion
- State-of-the-art infrared and thermal conductivity detectors with no moving parts and no manual adjustments
- User-friendly Cornerstone brand software
- Custom software interface designed specifically for touch-screen operation

LECO, Cornerstone are registered trademarks of LECO Corporation.





LECO Corporation 3000 Lakeview Avenue | St. Joseph, MI 49085 | 800-292-6141 | Phone: 269-985-5496 info@leco.com | www.leco.com | ISO-9001:2015 Q-994 | LECO is a registered trademark of LECO Corporation.

© 2019 LECO Corporation